表5. Xpert HBV VL分析的LLOQ的測定

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病毒	細菌	酵母菌	
BK Human polyoma virus	Human Immunodeficiency virus 1	<i>Staphylococcus epidermidis</i>	<i>Candida albicans</i>
Cytomegalovirus	Human Immunodeficiency virus 2	<i>Staphylococcus aureus</i>	
Epstein-Barr virus	Human papilloma virus 16		
Hepatitis A virus	Human papilloma virus 18		
Hepatitis C virus	Human T-cell lymphotropic virus type 1		
Herpes simplex virus 1	Human T-cell lymphotropic virus type 2		
Herpes simplex virus 2	Varicella Zoster virus		
Human herpes virus 6	Vaccina virus		
Human herpes virus 8			

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

- World Health Organization (WHO). Guidelines for the Prevention, Care and Treatment of Persons with Chronic Hepatitis B Infection. March 2015. Accessed March 14, 2018 at: <http://www.who.int/hiv/pub/hepatitis/hepatitis-b-guidelines/en/>.
- European Association for the Study of the Liver. EASL Clinical Practice Guidelines: Management of chronic hepatitis B virus infection. J. Hepat. 2012; 57:167-185. Available at: <http://dxdoi.org/10.1016/j.jhep.2012.02.010>.
- World Health Organisation. Global hepatitis report, 2017. WHO. April 2017.
- The 4th WHO International Standard for Hepatitis B Virus DNA for Nucleic Acid Amplification Techniques (NIBSC code: 10/266). National Institute for Biological Standards and Control; 2016.
- Centers for Disease Control and Prevention. Biosafety in Microbiological and Biomedical Laboratories (5th edition), accessed March 5, 2018 at <http://www.cdc.gov/biosafety/publications/>.
- Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections, Approved Guideline. Document M29 (refer to latest edition).
- World Health Organization. Safe management of wastes from health-care activities. 2nd Edition. WHO, 2014. Accessed April 20, 2018 at http://www.who.int/water_sanitation_health/publications/wastemanag/en/.
- REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on the classification labeling and packaging of substances and mixtures amending and repealing. List of Precautionary Statements, Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.
- Occupational Safety and Health Standards, Hazard Communication, Toxic and Hazard Substances (March 26, 2012) (29 C.F.R., pt. 1910, subpt. Z).
- Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Clinical and Laboratory Standards Institute, Wayne, PA, 2012.
- Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition. CLSI document EP05-A3. Clinical and Laboratory Standards Institute, Wayne, PA, 2014.
- Clinical and Laboratory Standards Institute. User Verification of Precision and Estimation of Bias; Approved Guideline – Third Edition. CLSI document EP15-A3. Clinical and Laboratory Standards Institute, Wayne, PA, 2014.
- Clinical and Laboratory Standards Institute. Evaluation of Linearity of Quantitative Measurement Procedures: A Statistical Approach. Approved Guideline. CLSI document EP06-A. Clinical and Laboratory Standards Institute, Wayne, PA, 2003.

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